Best Fit

08

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SISTEMA CLICK & GRIP

ISOGRIP®

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/ THE ISOGRIP® SYSTEM FOR **IMMEDIATE LOAD.**

COMPLETE STRUCTURE CREATION.

/ HOW IT WORKS

/STEP 01 FIXATION
Screw the **Provisional Expander Core (E)** to the Multi-Unit pillars with the **Fix Screw (B)**, and fit the **External Cylinders (A)** exerting slight pressure, until a “click” is felt.

/STEP 02 OVERFLOW
Overlay the prosthesis with resin.

/STEP 03 UNION
Union of the overlay prosthesis on the **Isogrip® Provisional System (E)** and solidification of the resin.

/STEP 04 SEPARATION
Separate the set -formed by the prosthesis including the **External Cylinders (A)** linked inside- and from the Multi-Unit pillars, exerting slight upward pressure.

/STEP 05 DRILLING
Fit the **Drilling Guides (H)** on the inside of the **External Cylinders (A)**, and insert the **Cutting Drill (G)** to drill the channels of the **Fixation Screws (B)** in this exact position.

/STEP 06 REMOVAL
Unscrew the **Provisional Expander Nucleus (E)**, from the Multi-Unit pillars.

/STEP 07 INSERTION
Insert one by one, the **Definitive Expansion Cores (D)** in the **External Cylinders (A)**, and make sure that the **Activating Screw (C)** does not exert pressure on the upper part of the **Definitive Expansion Core (D)**. Next, tighten the **Activation Screws (C)** with the **Activation Wrench (F)** in an counter clockwise.


/STEP 08 FINAL SCREWING
Position the **prosthesis** with the **Definitive Isogrip® System** included on the Multi-Unit pillars. Finish by screwing the complete structure using the **Fixing Screws (B)** to the implant.

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
L

/ OUTER CYLINDER

Bonds the Isogrip® System to the slipped prosthesis.

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Código: G09160

/ FIX SCREW

Responsible for joining the Isogrip® System to the Multi-Unit abutment.

Code: G09160

/ Carrying out immediate loads with the reversible Isogrip® System.

Carrying out immediate loads on any type of dental over Multi-Unit pilar, thanks to the Isogrip®.

This system allows:

- Make the screw channels in their exact position;
- Time saving;

In just 8 steps, we obtain a denture screwed on implants.

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/ ACTIVATION SCREW

Piece that houses the seat of the fixing screw; Activates the ultimate expanding core when tighten with the Activation Key (F).

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/ DEFINITIVE EXPANSION CORE

Piece that is attached to the provisional cylinder through an expansion of his upper area, produced by the pressure exerted by the activator screw.

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/ PROVISIONAL EXPANSION CORE

Piece geometrically analogous to the Isogrip abutment, with the only difference that it is machined in PEEK and does not require an activator screw.

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Código: E001564

/ ACTIVATION KEY

Key to tighten the **Activator Screw (C)**.

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Código: E000502

/ CUTTING DRILL

Drill required to drill the **Fix Screw (B) channels**.

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Código: E73001

/ DRILLING GUIDE

Guide to direct the cutter on this **exact position** and make the screw **channel super precisely**.

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Código: D30054

/ SCREWDRIVER

Used to screw the **Fix Screw (B)** to the **Isogrip® System**.

Code: D30054