



		Prime Ti-Base	Universal Ti-Base	CrCo-Base	Direct
Bego Semados®	Internal 3.25 / 3.75	X	X	✓	✓
Bego Semados®	Internal 4.1	X	X	✓	✓
Bego Semados®	Internal 4.5	X	X	✓	✓
Bego Semados®	Internal 5.5	X	X	✓	✓
Best Fit®	e-Clean Mini (055 HXE)	✓	✓	✓	✓
Best Fit®	e-Clean (134 HXE)	✓	✓	✓	✓
Best Fit®	e-Clean Wide (054 HXE)	X	✓	✓	✓
Best Fit®	s-Clean Narrow (010 HXI)	✓	X	✓	✓
Best Fit®	s-Clean (008 HXI)	✓	✓	✓	✓
Best Fit®	Internal Hex 002 HXI	✓	✓	✓	✓
Best Fit®	i-Clean (001 OCT)	✓	✓	✓	✓
Best Fit®	i-Clean over Octa Abutment	X	✓	✓	✓
Biohorizons®	Tapered 3.0	X	X	✓	✓
Biohorizons®	Tapered 3.5	✓	✓	✓	✓
Biohorizons®	Tapered 4.5	✓	✓	✓	✓
Biohorizons®	Tapered 5.7	X	X	✓	✓
Biomet 3i®	External 3.4 micro-mini®	X	✓	✓	✓
Biomet 3i®	External 4.1	✓	✓	✓	✓
Biomet 3i®	External 5.0	X	✓	✓	✓
Biomet 3i®	Certain® 3.4	X	✓	✓	✓
Biomet 3i®	Certain® 4.1	X	✓	✓	✓
Biomet 3i®	Certain® 5.0	X	✓	✓	✓
Bti®	External Tiny® 3.5	X	X	✓	✓
Bti®	External Universal® 4.1	✓	✓	✓	✓
Bti®	External Wide 5.5	X	✓	✓	✓
Bti®	Internal Narrow 3.5	X	X	✓	✓
Bti®	Internal Universal® 4.1	X	X	✓	✓
Bti®	Internal Wide 5.5	X	X	✓	✓
Camlog®	Camlog® System 3.3	X	X	✓	✓
Camlog®	Camlog® System 3.8	X	X	✓	✓
Camlog®	Camlog® System 4.3	X	X	✓	✓
Dentsply®	Astra® OsseoSpeed® 3.0 Yellow	✓	X	✓	✓
Dentsply®	Astra® OsseoSpeed® 3.5/4.0 Aqua	X	X	✓	✓
Dentsply®	Astra® OsseoSpeed® 4.5/5.0 Lilac	✓	✓	✓	✓
Dentsply®	Astra® Ossp. 20° Uni Abutment	✓	✓	✓	✓
Dentsply®	Astra® Ossp. 45° Uni Abutment	X	X	✓	✓
Dentsply®	Friadent® Int. Xive® 3.4	X	✓	✓	✓
Dentsply®	Friadent® Int. Xive® 3.8	X	✓	✓	✓
Dentsply®	Friadent® Int. Xive® 4.5	X	✓	✓	✓
Dentsply®	Astra® EV 33° Uni Abutment	X	X	✓	✓
Klockner®	SK2/NK2 4.2	X	✓	✓	✓
Klockner®	Essential®	X	✓	✓	✓
Klockner®	Octacone® Abutment	X	✓	✓	✓
Klockner®	Vega® MV 3.0	✓	X	✓	✓
Klockner®	Vega® NV 3.5	✓	X	✓	✓
Klockner®	Vega® RV 4.0	✓	X	✓	✓



		Prime Ti-Base	Universal Ti-Base	CrCo-Base	Direct
Mis®	Seven® Internal Narrow ø3.3	X	X	✓	✓
Mis®	Seven® Internal Standard ø3.5	✓	✓	✓	✓
Mis®	Seven® Internal Wide ø4.5	✓	✓	✓	✓
Mis®	Conical® C1/V3 Standard	X	X	✓	✓
Nobel Biocare®	Branemark® NP	✓	✓	✓	✓
Nobel Biocare®	Branemark® RP	✓	✓	✓	✓
Nobel Biocare®	Branemark® WP	X	✓	✓	✓
Nobel Biocare®	Active® NP	✓	✓	✓	✓
Nobel Biocare®	Active® RP	✓	✓	✓	✓
Nobel Biocare®	Replace® 3.5	X	✓	X	X
Nobel Biocare®	Replace® 4.3	X	✓	X	X
Nobel Biocare®	Replace® 5.0	X	✓	X	X
Nobel Biocare®	Multi-Unit® RP	✓	✓	✓	✓
Nobel Biocare®	Multi-Unit® WP	X	✓	✓	✓
Phibo®	External TSH® S2	X	✓	✓	✓
Phibo®	External TSH® S3/S4	X	X	✓	✓
Phibo®	External TSH® S5	X	✓	✓	✓
Phibo®	Aurea® Evo NP	X	X	✓	✓
Phibo®	Aurea® Evo RP	X	X	✓	✓
Phibo®	Aurea® Evo WP	X	X	✓	✓
Straumann®	Tissue Level RN®	✓	✓	✓	✓
Straumann®	Over Synocta® RN	X	✓	✓	✓
Straumann®	Tissue Level WN®	X	✓	✓	✓
Straumann®	Over Synocta® WN	X	✓	✓	✓
Straumann®	BLX® RB / WB	X	X	✓	✓
Straumann®	Bone Level NC®	✓	✓	✓	✓
Straumann®	Bone Level RC®	✓	✓	✓	✓
Straumann®	Over Screw Retained Abutment	X	X	✓	✓
Sweden & Martina®	External Outlink2® ø3.75/4.10 y 5.0	✓	✓	✓	✓
Sweden & Martina®	Int. Premium® 3.3	X	X	✓	✓
Sweden & Martina®	Int. Premium® y Kohno® 3.8	X	X	✓	✓
Sweden & Martina®	Int. Premium® y Kohno® 4.25	X	X	✓	✓
Sweden & Martina®	Int. Premium® y Kohno® 5.0	X	X	✓	✓
Zimmer®	Tapered Screw Vent® 3.5	✓	✓	✓	✓
Zimmer®	Tapered Screw Vent® 4.5	✓	✓	✓	✓
Zimmer®	Tapered Screw Vent® 5.7	X	✓	✓	✓